

icoone®

BEAUTY

Protocol Cards




Micro Actions, Macro Changes

INTRODUCTION

This document is intended to accompany all operators that use icoone around the world.
It will provide protocols to be used with your customers, in particular for the creation of treatment courses.

However, we would like to remind you that, while following the guidelines provided, each protocol must be carefully customized, after visual and palpatory evaluation of the client's silhouette and based on his specific requests.

We are sure that our technology, combined with your professionalism, will allow you to achieve exceptional results. Good work!



PROTOCOLS

Perfect shape
Immune symmetry
Silky appearance
Firm solution
Mama care
Time stop

For Women & Men
Body, Face, Neck
and Décolleté

NUMBER OF SESSIONS

5 sessions
10 sessions or more

The results are
evident even after
the first
5 sessions

TYPE OF PACKAGES

EXPRESS

Sessions every day
for faster results

BASE

2-3 sessions per
week

MAINTENANCE

1-2 sessions per
month after obtaining
the desired result

STRUCTURAL LOGIC OF BEAUTY PROTOCOLS

1° STEP : GENERAL TOTAL BODY PREPARATION

General homogeneous skin condition
all over the body:

1 basic program (20')



20"

Two different skin conditions
(e.g. above and below):

2 different basic programs (10' + 10')



10"



10"

2° STEP: LOCALIZED TREATMENT

Many focus programs (10 min each), depending on the localized aesthetic problems
you intend to treat

Example: Abdoton (10') + Modelplus (10')
If the focus of the treatment is on toning
the abdomen and obtaining a push-up
effect at the level of the buttocks.



10"



10"

BEAUTY PROTOCOLS FOR BODY, FACE, NECK AND DÉCOLLETÉ



REMODELING:

Effective on the most resistant adipose deposits, it is particularly suitable for targeting localized imperfections and specific areas such as hips, abdomen, arms, legs, knees, calves, back, and buttocks.



ANTI - CELLULITE:

It contributes to smoothing the skin, improving fluid circulation, and harmonizing body shapes.



FIRMING AND TONING:

Promotes the stimulation of collagen and elastin, restoring the natural firmness of tissues. Particularly recommended in cases of skin laxity due to weight loss.



LYMPHATIC DRAINAGE:

Microstimulation allows acting globally on the body and face, reactivating circulation and the natural functioning of the body. Thanks to lymphatic drainage, tissues appear visibly regenerated, and the skin revitalized.



ANTI-AGING:

Effectively targets signs of aging such as wrinkles and sagging skin. It is such a gentle and non-invasive treatment that can be done daily for quicker and long-lasting results.



PREGNANCY:

Offers expectant mothers the opportunity to experience an effective and relaxing treatment. By draining fluids, it helps to relieve leg discomfort, control water retention, and prevent the appearance of imperfections such as cellulite and stretch marks.

INDEX

1. Perfect shape	
1.1. Full body	1
1.2. Waist line	3
1.3. Upper body	5
1.4. Lower body	7
1.5. Face	9
2. Firm Solution	
2.1. Full body.....	11
2.2. Arms	13
2.3. Abdomen	15
2.4. Inner thigh.....	17
2.5. Buttocks	19
2.6. Specific action	21
2.8. Face	23
3. Time Stop	
3.1. Face	25
4. Immune Symmetry	
4.1.Venous drainage	27
4.2. Lymphatic drainage	28
4.3. Face drainage	29
5. Silky Appearance	
5.1. Fibrotic Cellulite.....	30
5.2. Adipose Cellulite	32
5.3. Edematous Cellulite	34
6. Mamma care	
6.1. Mamma care	36

FULL BODY 5

Designed to treat the whole body with a specific remodeling action on particularly difficult and resistant fatty areas. A deep lipolytic action accompanied by a pleasant feeling of well-being and lightness.

TARGET

Women
Men*

FREQUENCY

Twice a week

SESSION DURATION

60 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Dry Skin
Thick Skin
Skin with adipose accumulations

BODY PART

Abdomen, waist, hips, thighs, legs, ankles, arms, back, hump of a bison

PROTOCOL

FROM SESSION 1 - 5

BASIC



CELLFAT 20"



FOCUS



SILHOUETTE 10" MENZONE 10" HIGHPLUS 10" REMOD 10"

FULL BODY 10

Designed to treat the whole body with a specific remodeling action on particularly difficult fatty areas. An even more effective and profound lipolytic action to treat resistant and advanced adipose overload.

TARGET

Women
Men*

FREQUENCY

Twice a week

SESSION DURATION

60 - 70 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Dry Skin
Thick Skin
Skin with adipose accumulations

BODY PART

Abdomen, waist, hips, thighs, legs, ankles, arms, back, hump of a bison

PROTOCOL

FROM SESSION 1 - 5

BASIC



CELLFAT 20"



FOCUS



SILHOUETTE 10" MENZONE 10" HIGHPLUS 10" REMOD 10"

FROM SESSION 6 - 10

POWER UP



CELLDRAIN 20"



SILHOUETTE 10" MENZONE 10" HIGHPLUS 10" REMOD 10" FATZONE 10"

WAIST LINE 5

Designed to treat adipose accumulations in the waist area before they spread. A deep lipolytic action accompanied by a pleasant feeling of well-being and lightness.

TARGET

Women
Men*

FREQUENCY

Twice a week

SESSION DURATION

40 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Dry Skin
Thick Skin
Skin with adipose
accumulations

BODY PART

Abdomen
Waist

PROTOCOL

FROM SESSION 1 - 5

BASIC



CELLFAT 20"



FOCUS



SILHOUETTE 10"



MENZONE 10"

WAIST LINE 10

Designed to treat significant adipose accumulations in the waist area both in terms of volume and extension of the area. A deep lipolytic action accompanied by a pleasant feeling of well-being and lightness.

TARGET

Women
Men*

FREQUENCY

Twice a week

SESSION DURATION

50 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Dry Skin
Thick Skin
Skin with adipose accumulations

BODY PART

Abdomen
Waist

PROTOCOL

FROM SESSION 1 - 5

BASIC



CELLFAT 20"



FOCUS



SILHOUETTE 10"



MENZONE 10"

FROM SESSION 6 - 10

POWER UP



CELLDRAIN 20"



SILHOUETTE 10"



MENZONE 10"



FATZONE 10"

UPPER BODY 5

Designed to treat adipose accumulations in the upper part of the body, specifically in the arm area, before they spread. A deep lipolytic action accompanied by a pleasant feeling of well-being and lightness.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Dry Skin
Thick Skin
Skin with adipose accumulations

BODY PART

Arms
Schiena
Bison hump

PROTOCOL

FROM SESSION 1 - 5

BASIC



CELLFAT 20"



FOCUS



HIGHPLUS 10"

UPPER BODY 10

Designed to treat in the upper part of the body, specifically in the arm area, important adipose accumulations both in terms of volume and extension of the area. A deep lipolytic action accompanied by a pleasant feeling of well-being and lightness.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Dry Skin
Thick Skin
Skin with adipose accumulations

BODY PART

Arms
Schiena
Bison hump

PROTOCOL

FROM SESSION 1 - 5

BASIC



CELLFAT 20"



FOCUS



HIGHPLUS 10"

FROM SESSION 6 - 10

POWER UP



CELLFAT 20"



ARMTON 10"



HIGHPLUS 10"

LOWER BODY 5

Designed to treat adipose accumulations in the lower part of the body, from the waist down, before they spread. A deep lipolytic action accompanied by a pleasant feeling of well-being and lightness.

TARGET

Women
Men*

FREQUENCY

Twice a week

SESSION DURATION

40 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Dry Skin
Thick Skin
Skin with adipose accumulations

BODY PART

Hips
Thighs
Legs
Ankles

PROTOCOL

FROM SESSION 1 - 5

BASIC



CELLFAT 20"



FOCUS



SILHOUETTE 10"



REMOD 10"

LOWER BODY 10

Designed to treat important adipose accumulations in the lower part of the body, from the waist down, both in terms of volume and extension of the area. A deep lipolytic action accompanied by a pleasant feeling of well-being and lightness.

TARGET

Women
Men*

FREQUENCY

Twice a week

SESSION DURATION

50 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Dry Skin
Thick Skin
Skin with adipose accumulations

BODY PART

Hips
Thighs
Legs
Ankles

PROTOCOL

FROM SESSION 1 - 5

BASIC



CELLFAT 20"



FOCUS



SILHOUETTE 10"



REMOD 10"

FROM SESSION 6 - 10

POWER UP



CELLDRAIN 20"



SILHOUETTE 10"



CELLFAT** 10"



REMOD 10"

FACE 5

Designed to treat the double chin before adipose accumulations spread.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Dry Skin
Thick Skin
Skin with adipose
accumulations

BODY PART

Double chin

PROTOCOL

FROM SESSION 1 - 5

BASIC



NECK FACE DRAIN 20"



FOCUS



CHIN SCULPT 10"

FACE 10

Designed to treat the chin with significant adipose accumulations both in terms of volume and extension of the area.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

40 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Dry Skin
Thick Skin
Skin with adipose
accumulations

BODY PART

Double chin

PROTOCOL

FROM SESSION 1 - 5

BASIC



NECK FACE DRAIN 20"



FOCUS



CHINSCULPT 10"

FROM SESSION 6 - 10

POWER UP



NECK FACE DRAIN 20"



NECK FACE ONE 10"



CHINSCULPT 10"

FULL BODY 5

Ideal for an immediate firming effect on the skin all over the body. Thanks to the effective action of stimulating the production of collagen and elastin and the use of ergonomic handpieces, it is possible to specifically treat difficult and often delicate areas such as the inner thighs, the inner arms, the abdomen and the buttocks.

TARGET

Women
Men*

FREQUENCY

Twice a week

SESSION DURATION

60 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Whole body

PROTOCOL

FROM SESSION 1 - 5

BASIC



BIOYOUNG 20"



FOCUS



ABDOTON 10"



ARMTON 10"



INTON 10"



MODELPLUS 10"

FULL BODY 10

Ideal for an immediate firming effect on particularly loose and thin body skin. Thanks to the effective action of stimulating the production of collagen and elastin and the use of ergonomic handpieces, it is possible to specifically treat difficult and often delicate areas such as the inner thighs, the inner arms, the abdomen and the buttocks.

TARGET

Women
Men*

FREQUENCY

Twice a week

SESSION DURATION

60 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Whole body

PROTOCOL

FROM SESSION 1 - 5

BASIC



BIOYOUNG 20"



FOCUS



ABDOTON 10"



ARMTON 10"



INTON 10"



MODELPLUS 10"

FROM SESSION 6 - 10

POWER UP



ELASTO 20"



ABDOTON 10"



ARMTON 10"



INTON 10"



MODELPLUS 10"

*NO INTON AND MODELPLUS PROGRAMS, INSERT FOCUS, BIOYOUNG AS A THIRD FOR 10 MIN IN AREAS WHERE NECESSARY.

THE PROTOCOLS DESCRIBED WERE STRUCTURED ON THE BASIS OF VISUAL AND PALPATORY EVALUATION. EACH PROTOCOL CAN BE ADAPTED AND CUSTOMIZED ACCORDING TO THE CUSTOMER'S NEEDS

ARMS 5

Ideal for an immediate firming effect on the skin of the whole body that gives at the same time a pleasant feeling of well-being and lightness. Thanks to the effective action of stimulating the production of collagen and elastin and the use of ergonomic handpieces, it is possible to specifically treat difficult areas such as the arms.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Arms

PROTOCOL

FROM SESSION 1 - 5

BASIC



BIOYOUNG 20"



FOCUS



ARMTON 10"

ARMS 10

Ideal for an immediate firming effect on the whole body in the presence of particularly loose and delicate skin. Thanks to the effective action of stimulating the production of collagen and elastin and the use of ergonomic handpieces, it is possible to specifically treat difficult areas such as the arms.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

40 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Arms

PROTOCOL

FROM SESSION 1 - 5

BASIC



BIOYOUNG 20"



FOCUS



ARMTON 10"

FROM SESSION 6 - 10

POWER UP



BIOYOUNG 20"



ARMTON 10"



ELASTO* 10"

ABDOMEN 5

Ideal for an immediate firming effect on the skin all over the body. Thanks to the effective action of stimulating the production of collagen and elastin and the use of ergonomic handpieces, it is possible to specifically treat difficult areas such as the abdomen.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

40 minutes

TREATMENT SESSIONS

5 sedute

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Arms

PROTOCOL

FROM SESSION 1 - 5

BASIC



BIOYOUNG 20"



FOCUS



ABDOTON 10"

ABDOMEN 10

Ideal for an immediate firming effect on the skin all over the body. Thanks to the effective action of stimulating the production of collagen and elastin and the use of ergonomic handpieces, it is possible to specifically treat difficult areas such as the abdomen even after significant weight loss.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Arms

PROTOCOL

FROM SESSION 1 - 5

BASIC



BIOYOUNG 20"



FOCUS



ABDOTON 10"

FROM SESSION 6 - 10

POWER UP



ELASTO 20"



ABDOTON 10"

INNER THIGH 5

Ideal for an immediate firming effect on the skin all over the body. Thanks to the effective action of stimulating the production of collagen and elastin and the use of ergonomic handpieces, it is possible to specifically treat difficult areas such as the inner thighs.

TARGET

Women

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Inner thigh

PROTOCOL

FROM SESSION 1 - 5

BASIC



BIOYOUNG 20"



FOCUS



INTON 10"

INNER THIGH 10

Ideal for an immediate firming effect on particularly loose and delicate skin all over the body. Thanks to the effective action of stimulating the production of collagen and elastin and the use of ergonomic handpieces, it is possible to specifically treat difficult areas such as the inner thighs.

TARGET

Women

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Inner thigh

PROTOCOL

FROM SESSION 1 - 5

BASIC



BIOYOUNG 20"



FOCUS



INTON 10"

FROM SESSION 6 - 10

POWER UP



ELASTO 20"



INTON 10"

BUTTOCKS 5

It firms the skin of the body by virtue of the effective action of stimulating the production of collagen and elastin. Thanks to the use of ergonomic handpieces, it is possible to treat all parts of the body, even the most difficult ones. Ideal for a tonic and anti-gravity effect on the buttocks.

TARGET

Women

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Buttocks

PROTOCOL

FROM SESSION 1 - 5

BASIC



BIOYOUNG 20"



FOCUS



MODELPLUS 10"

BUTTOCKS 10

It firms the skin of the body, particularly loose and thin, by virtue of the effective action of stimulating the production of collagen and elastin. Thanks to the use of ergonomic handpieces, it is possible to treat all parts of the body, even the most difficult ones. Ideal for a tonic and anti-gravity effect on the buttocks.

TARGET

Women

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Buttocks

PROTOCOL

FROM SESSION 1 - 5

BASIC



BIOYOUNG 20"



FOCUS



MODELPLUS 10"

FROM SESSION 6 - 10

POWER UP



ELASTO 20"



MODELPLUS 10"

SPECIFIC ACTION 5

Ideal for an immediate firming effect on the skin thanks to the effective action of stimulating the production of collagen and elastin. By virtue of the use of ergonomic handpieces, it is possible to treat all parts of the body, those that are smaller, more difficult and delicate, such as the inner thighs, knees and breasts, and those that have stretch marks.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Breast
Kneecchia
Stretch marks

PROTOCOL

FROM SESSION 1 - 5

BASIC



ELASTO 20"



FOCUS



BIOYOUNG 10"

SPECIFIC ACTION 10

Ideal for an immediate firming effect on the skin thanks to the effective action of stimulating the production of collagen and elastin. By virtue of the use of ergonomic handpieces, it is possible to treat all parts of the body, those that are smaller, more difficult and delicate, such as the inner thighs, knees and breasts, and those that have stretch marks.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Breast
Kneecchia
Stretch marks

PROTOCOL

FROM SESSION 1 - 5

BASIC



ELASTO 20"



FOCUS



BIOYOUNG 10"

FROM SESSION 6 - 10

POWER UP



ELASTO 20"



BIOYOUNG 10"

FACE 5

Designed to treat the face in the presence of relaxed skin for an immediate effect of general turgor and luminosity.

TARGET

Women

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Face
Neck

PROTOCOL

FROM SESSION 1 - 5

BASIC



NECK FACE DRAIN 20"



FOCUS



NECK FACE ONE 10"

FACE 10

Designed to treat the face in the presence of very relaxed skin for an immediate lifting and firming effect.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Face
Neck

PROTOCOL

BASIC



FOCUS

FROM SESSION 1 - 5



NECK FACE DRAIN 20"



NECK FACE ONE 10"

FROM SESSION 6 - 10

POWER UP



NECK FACE ONE 20"



NECK FACE TWO 10"

FACE 5

Ideal for treating the face in the presence of relaxed skin and the presence of wrinkles in the orbicular area of the mouth and eyes.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Face
Neck

PROTOCOL

FROM SESSION 1 - 5

BASIC



FOCUS



NECK FACE ONE 20"



SMOOTH FACE 10"

FACE 10

Ideal for treating the face in the presence of very relaxed and thin skin, with marked wrinkles in the orbicular area of the mouth and eyes.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Flabby Skin
Thin skin
Skin with loss of tone

BODY PART

Face
Neck

PROTOCOL

BASIC



FOCUS

FROM SESSION 1 - 5



NECK FACE ONE 20"



SMOOTH FACE 10"

FROM SESSION 6 - 10

POWER UP



NECK FACE ONE 10"



NECK FACE TWO 10"



SMOOTH FACE 10"

VENOUS DRAINAGE

Designed to support the proper functioning of the immune system thanks to the action of stimulating the capillary network and the consequent improvement of venous circulation. It promotes comfort and lightness to the legs and redefines the ankles. It can be used to restore good health after a period of stress.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

40 minutes

TREATMENT SESSIONS

5 sessions
10 sessions

PHYSICAL EXAMINATION

Dry skin
Dehydrated Skin
Skin with weak capillaries

BODY PART

Full body

BASIC



FOCUS

PROTOCOL

5 SESSIONS



FLOWING 20"



CELLDRAIN* 10"



REM0D 10"

10 SESSIONS

POWER UP



FLOWING 20"



CELLDRAIN* 10"



REM0D 10"

LYMPHATIC DRAINAGE

Designed to support the proper functioning of the immune system thanks to the action of stimulating lymphatic drainage with the consequent elimination of water retention. It can be used for preventive purposes to maintain a state of well-being.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

50 minutes

TREATMENT SESSIONS

5 sessions
10 sessions

PHYSICAL EXAMINATION

Dry skin
Dehydrated Skin
Skin with weak capillaries

BODY PART

Full body

BASIC



FOCUS

5 SESSIONS



CELLDRAIN 20''



REMOD 10''



TONUS* 10''



OPTIMUM* 10''

10 SESSIONS

POWER UP



CELLDRAIN 20''



REMOD 10''



TONUS* 10''



OPTIMUM* 10''

PROTOCOL

FACE DRAINAGE

Designed to promote fluid drainage and give the face a state of new freshness and luminosity, deflating it, detoxifying the skin and reducing bags and dark circles under the eyes.

TARGET

Women
Men

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

5 sessions
10 sessions

PHYSICAL EXAMINATION

Dry skin
Dehydrated Skin
Skin with weak capillaries

BODY PART

Face
Neck
Décolleté

BASIC



FOCUS

5 SESSIONS



NECK FACE DRAIN 20"



EYELIGHT 10"

10 SESSIONS

POWER UP



NECK FACE DRAIN 20"



EYELIGHT 10"

PROTOCOL

FIBROTIC CELLULITE 5

Ideal for softening accumulations and improving fluid circulation, thus reducing the appearance of recently formed fibrous cellulite.

TARGET

Women

FREQUENCY

Twice a week

SESSION DURATION

40 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Thick Skin
Dehydrated Skin
Dry Skin
Orange Peel Skin

BODY PART

Lower body

PROTOCOL

FROM SESSION 1 - 5

BASIC



SILK 20"



FOCUS



CELLDRAIN* 10"



SILHOUETTE 10"

FIBROTIC CELLULITE 10

Ideal for softening accumulations and improving fluid circulation, thus reducing the appearance of old orange peel skin. with the consequent improvement of fluid circulation and the restoration of harmonious shapes.

TARGET

Women

FREQUENCY

Twice a week

SESSION DURATION

50 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Thick Skin
Dehydrated Skin
Dry Skin
Orange Peel Skin

BODY PART

Lower body

PROTOCOL

FROM SESSION 1 - 5

BASIC



SILK 20"



FOCUS



CELLDRAIN* 10"



SILHOUETTE 10"

FROM SESSION 6 - 10

POWER UP



ELASTO 20"



SILK* 10"



MODELPLUS 10"



FATZONE 10"

*USE ON LEGS 10MIN

ADIPOSE CELLULITE 5

Ideal for treating adipose accumulations mixed with cellulite, especially in the waist, buttocks and thighs. It improves the circulation of liquids, reduces the volume of fat with the consequent restoration of harmonious shapes.

TARGET

Women

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Thick Skin
Dehydrated Skin
Dry Skin
Skin with adipose accumulations

BODY PART

Lower body

PROTOCOL

FROM SESSION 1 - 5

BASIC



CELLFAT 20"



FOCUS



SILHOUETTE 10"

A DIPOSE CELLULITE 10

Ideal for treating long-standing adipose accumulations mixed with cellulite. It improves the circulation of liquids, reduces the volume of fat with the consequent restoration of harmonious shapes.

TARGET

Women

FREQUENCY

Twice a week

SESSION DURATION

50 minutes

TREATMENT SESSIONS

10 sessions (5+5)

PHYSICAL EXAMINATION

Thick Skin
Dehydrated Skin
Dry Skin
Skin with adipose accumulations

BODY PART

Lower body

PROTOCOL

FROM SESSION 1 - 5

BASIC



CELLFAT 20"



FOCUS



SILHOUETTE 10"

FROM SESSION 6 - 10

POWER UP



CELLDRAIN 20"



CELLFAT* 10"



REMOD 10"



FATZONE 10"

*USE ON LEGS 10MIN

EDEMATOUS CELLULITE 5

Ideal for improving circulation, thus reducing newly formed water retention, with the consequent restoration of a harmonious appearance.

TARGET

Women

FREQUENCY

Twice a week

SESSION DURATION

30 minutes

TREATMENT SESSIONS

5 sessions

PHYSICAL EXAMINATION

Dehydrated Skin
Dry Skin
Edema due to the
presence of cellulite

BODY PART

Lower body

PROTOCOL

FROM SESSION 1 - 5

BASIC



FLOWING 20"



FOCUS



CELLDRAIN* 10"

EDEMATOUS CELLULITE 10

Ideal for improving circulation, thus reducing long-standing water retention located in the waist and thighs. The result is the restoration of a harmonious appearance.

TARGET

Women

FREQUENCY

Twice a week

SESSION DURATION

40 minutes

SEDUTE CONSIGLIATA

10 sessions (5+5)

ESAME FISICO

Dehydrated Skin
Dry Skin
Edema due to the
presence of cellulite

AREE DI RIFERIMENTO

Lower body

PROTOCOL

FROM SESSION 1 - 5

BASIC



FLOWING 10"



FOCUS



CELLDRAIN* 10"

FROM SESSION 6 - 10

POWER UP



CELLDRAIN 20"



FLOWING* 10"



SILHOUETTE 10"

MAMA CARE

For all women who want to prevent or minimize pregnancy** related imperfections, to maintain the harmony of body shapes and good skin quality in the affected areas.

TARGET

Women

FREQUENCY

Twice a week

SESSION DURATION

50 minutes

SEDUTE CONSIGLIATA

5 sessions
10 sessions

ESAME FISICO

Dehydrated Skin
Dry Skin
Skin with weak
capillaries

AREE DI RIFERIMENTO

Total body
(Abdomen, lower back)

BASIC



FOCUS

PROTOCOL

5 SESSIONS



FLOWING 10"



CELLDRAIN 10"



SILK* 10"



TONUS* 10"



REMOD 10"

10 SESSIONS

POWER UP



FLOWING 10"



CELLDRAIN 10"



SILK* 10"



TONUS* 10"



REMOD 10"